



CHARLOTTESVILLE-ALBEMARLE
RESCUE SQUAD



Glucose Monitor Weekly Control Documentary

(3 weeks per sheet)

Bayer Ascensia Contour Blood Glucose Monitoring System

Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
<i>Ex: 10/9/07</i>	<i>Ex: #1</i>	<i>Ex: 3078148</i>	<i>Ex: 7AC3C53</i>	<i>Ex: 01/09</i>	<i>Ex: 100-138</i>
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service
<i>Ex: 1508</i>	<i>Ex: 03/09</i>	<i>Ex: 107</i>	<i>Ex: Pass</i>	<i>Ex: N/A</i>	<i>Ex: N/A</i>
XXXXX	XXXXX	XXXXX	XXXXX	XXXXX	XXXXX
Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service
XXXXX	XXXXX	XXXXX	XXXXX	XXXXX	XXXXX
Date of Test	CARS Device #	BGMS Serial #	Test Strip Lot #	Strip Expiration Date	Normal Control Range
Normal Control Solution Lot #	Control Solution Expiration Date	Test Result	Pass or Fail? If Fail Repeat.	Repeated Test Result	Fail 2nd time? YES, Take Out of Service